Please read this package insert carefully before use.

English

Capilia™ TB-Neo

INTENDED USE

To detect *M. tuberculosis* complex antigens in a suspension of bacteria cultured on the medium for acid-fast bacteria (AFB) or in a portion of the liquid medium for AFB in which bacteria have been cultured (to assist in the diagnosis of tuberculosis).

SUMMARY AND EXPLANATION OF THE TEST

Tuberculosis ("TB") is a chronic disease caused by infection with *Mycobacterium tuberculosis*. The annual number of new TB cases is around 8.8 million, and the number of deaths from TB is around 1.5 million.

The number of infections by non-tuberculous mycobacteria ("NTM") has been increasing in recent years. NTM symptoms can be similar to those caused by TB. It is critical to distinguish NTM infections from TB to determine the appropriate treatment because many NTM are resistant to ordinary antituberculous drugs.

Definite diagnosis of TB is made by detecting *M. tuberculosis* from a clinical sample.

The results for the conventional detection method of the *M. tuberculosis* complex (including the niacin accumulation test) usually require 4 to 8 weeks.


Recently, a nucleic acid probe-based identification method has been developed as a tool for rapid detection of the *M. tuberculosis* complex. However, it requires cumbersome handling procedures, special instruments/equipment and skilled operators.

Capilia TB-Neo adopts an immunochromatography method, which can detect the MPB64 antigens specifically produced by the *M. tuberculosis* complex.

Capilia TB-Neo is able to detect the *M. tuberculosis* complex in bacterial isolates specifically and with a high degree of sensitivity, and to provide rapid test results by a simple operation. No special instruments or equipment are required.

PRINCIPLE OF THE TEST

This product comprises a test plate with a carrier strip containing a sample placement area, a reagent area including a colloidal gold labeled anti-MPB64 monoclonal antibody (mouse) (hereafter referred to as "colloidal gold labeled anti-MPB64 antibody") and a developing area that fixes the anti-MPB64 monoclonal antibody (mouse) (hereafter referred to as "anti-MPB64 antibody").

When a sample is placed on the sample placement area of the test plate, the colloidal gold labeled anti-MPB64 antibody dissolves and forms an immune complex with MPB64 antigens in the sample. This immune complex migrates through the developing area by capillary action, is captured by the anti-MPB64 antibody fixed in the developing area, and forms a red purple line of colloidal gold in the reading area [T].

The red purple line visually displays the existence of MPB64 antigens in the sample.

Regardless of the existence of the MPB64 antigens in the sample, excess colloidal gold labeled anti-MPB64 antibody further migrates through the developing area, are captured by anti-mouse immunoglobulin antibodies fixed in the developing area, and form a red purple line in the reading area [C]. This means the colloidal gold labeled anti-MPB64 antibodies have migrated normally.

REAGENTS AND MATERIALS PROVIDED

[REF] CATB0870 Capilia TB-Neo (100 Tests)

[REF] CATB0871 Capilia TB-Neo (10 Tests)

Test plates

• Components

Colloidal gold labeled anti-MPB64 monoclonal antibody (mouse)
Anti-MPB64 monoclonal antibody (mouse)

MATERIALS REQUIRED BUT SOLD SEPARATELY

[REF] CATB0877 Capilia TB-Neo Extraction Buffer (20 mL/Bottle)

MATERIALS REQUIRED BUT NOT PROVIDED

Timer, micropipette, pipette tips

WARNING AND PRECAUTIONS

1. Precautions when handling samples

- 1) Use a suspension of bacteria cultured on the medium for AFB or a bacterial culture liquid as a sample. No clinical samples such as a human body fluid, tissue or bronchial lavage fluid can be directly used as they are as samples.
- 2) Cultured samples must be used immediately for testing.^{Note}
- 3) Bacteria inoculated in the media for AFB must be carefully handled as they may cause infections.
- 4) Some samples including substrains of *Mycobacterium bovis* BCG among the *M. tuberculosis* complex (Copenhagen, Glaxo, Pasteur, and Tice) are interpreted as negative, because no MPB64 is produced from these substrains.
- 5) Although no cross reactivity with 4 strains of *M. marinum*, which are described in the "Cross reactivity" section, was observed, cross reactivity with other strains of *M. marinum* has not yet been studied.
- 6) This product is unable to differentiate each of the *M. tuberculosis*, *M. bovis*, *M. africanum* and *M. microti*, from the others even if the *M. tuberculosis* complex test results are positive.
- 7) Strains with a mutated MPB64 gene may not be detected by this product.

^{Note} Positive cultures in liquid media and colonies on solid media can be tested for as long as one year when stored at -80°C for testing.

2. Precautions when using

- 1) This product is a rapid test for identification of *M. tuberculosis* complex.
Any clinical diagnosis should be made by an attending physician, in combination with the clinical symptoms, and other test results.
- 2) This product should be used in accordance with the procedure stated in this package insert.
- 3) In order to prevent deterioration, this product should be stored between 2°C and 30°C, avoiding high temperatures, high humidity and direct sunlight.

- 4) If this product has been refrigerated, it must be removed from the refrigerator at least 30 minutes before use to be acclimatized to room temperature.

- 5) The aluminum pouches containing the test plates should not be opened until they are about to be used.

- 6) The sample placement area and the reading area of the test plate should not be touched with the hands.

3. Precautions for disposal

- 1) Because used test plates, pipette tips, remaining samples and other test materials/devices used for identification of the *M. tuberculosis* complex may cause infections, they should be autoclaved (at 121°C, for over 30 min) before disposal.
- 2) When sterilized used reagents or other test materials are disposed of, they should be treated in accordance with the laws and regulations concerning medical waste disposal and water pollution control.

STORAGE CONDITIONS

Storage : Store at 2°C to 30°C.

Keep away from direct sunlight.

Do not use test plate beyond the expiration date.

DO NOT FREEZE.

SAMPLE COLLECTION AND PREPARATION

1. Methods of sample collection

- 1) Use a suspension of bacteria cultured on the medium for AFB or a bacterial culture liquid as a sample. No clinical samples such as a human body fluid, tissue or bronchial lavage fluid can be directly used, as they are as samples.
- 2) The testing laboratory that handles *M. tuberculosis* must be isolated from the external space by dual doors and the mycobacteria must be handled in a biological safety cabinet.
- 3) Bacteria inoculated into the media for AFB must be carefully handled, as they may cause infections.
- 4) Dispense samples with a micropipette and use a new tip with a filter for dispensing each sample.

2. Sample preparation

After completing appropriate pretreatment^{Note} of clinical samples such as human body fluids, tissue and bronchial lavage fluid, inoculate the sample into a culture medium for AFB.

^{Note} Example of pretreatment of sputum:

• N-acetyl-L-cysteine-sodium hydroxide (NALC-NaOH) method:

Add at least a 2-fold volume of NALC-NaOH solution to the container containing the collected sputum, stir the container with a vortex mixer and invert the container to expose the inside of the screw cap and the inner wall of the tube to the NALC-NaOH solution. Allow the container to stand at room temperature for 15 min and intermittently shake it lightly by hand. After adding a cold-sterilized, 10-fold volume of phosphate-buffered (PB) solution (pH 6.8) to the container and mixing it, centrifuge 3,000 × g for 20 min. Suspend the precipitate in 1 mL of PB solution. Then, inoculate 0.1 mL of the suspension into a 1% Ogawa medium or inoculate 0.5 mL of the same into a liquid medium.

• 4% sodium hydroxide method:

Add at least a 2-fold volume of 4% sodium hydroxide solution to the container containing the collected sputum, fully homogenize the mixture and immediately inoculate 0.1 mL of the mixture into a 3% Ogawa medium.

- 1) In the case of using a **liquid medium** for AFB (Example: Middle Brook 7H9 broth)
Incubate at 37°C for 1 to 3 weeks until the liquid medium becomes cloudy, because of the growth of bacteria. In the event that a Mycobacteria Growth Indicator Tube (MGIT) is used, incubate until a positive interpretation is possible. In both cases, it is necessary to confirm the presence of AFB by acid-fast staining. Stir the liquid medium in the culture tube and use the medium as a sample.
- 2) In the case of using a **solid medium** for AFB (Example: Ogawa medium)
Incubate at 37°C for 2 to 4 weeks until the growth of bacterial colonies is confirmed on the solid medium, and then confirm the presence of AFB by acid-fast staining.
 - ① Dispense 0.2mL of extraction buffer (sold separately) into a tube.
 - ② Remove 1 μL (equivalent to the amount of a 1 mm-diameter platinum micro-loop) from the bacterial colony that has grown on the solid medium.
 - ③ Suspend the collected bacteria in the buffer solution in the tube.
 - ④ Close the tube with a stopper and fully suspend with a vortex mixer. Then, use the bacterial suspension for samples.

TEST PROCEDURE

- 1) Dispense 80–100 μL of sample onto the sample placement area of the test plate.
- 2) Observe the reading area of the test plate after 15 min and interpret the result according to the "READING TEST RESULTS."

READING THE RESULTS

Allow the samples to react according to the procedure and read the red purple lines that appear in the reading area.

Positive		When red purple lines are seen at both [T] and [C] in the reading area (two lines), the result is read as positive. When a very faint red purple line is seen at [T] in the reading area, the result is also interpreted as positive.
Negative		When no red purple line is seen at [T] in the reading area but a red purple line is seen only at [C] in the reading area (one line), the result is read as negative. When a red purple line at [C] in the reading area is faint but visually recognizable, chromatographic development has occurred normally.
Retesting		When no red purple line is seen at [C] in the reading area, there may be some problem with the test procedure or the reagent quality. The test should be performed again, using another test plate.

If the line falls within any of the sections of the reading area, at whatever point, the test is considered valid. (The boundaries between the sections of the reading area are shown by the nicks in the reading window frame.)

(Note)

1. If no line is seen at [C] in the reading area, there may be some problem with the test procedure or the reagent quality. Therefore, the test should be performed again, using another test plate.
2. Do not use the test plate for a reading result beyond the judgement time (15 minutes) as result may change due to drying, etc.
3. If a test result using this product is interpreted as positive, the presence of the *M. tuberculosis* complex in the sample is strongly suggested. However, there is a possibility of combined infections of *M. tuberculosis* and NTM, and an infection with NTM cannot be excluded.
4. The sample placement area or the reading area of the test plate should not be touched with the hands.
5. Dispense samples with a micropipette and use a new tip with a filter for dispensing each sample.

6. Even if a test result using this product is interpreted as negative, this may be due to its inability to detect the *M. tuberculosis* complex when the MPB64 concentration in the sample is below the detection limit or a mutation arises in the MPB64 gene of the *M. tuberculosis* complex. ⁽²²⁾ Therefore, a negative result does not necessarily rule out the possibility of infection with *M. tuberculosis*.

Any clinical diagnosis must be made by an attending physician, in combination with the clinical symptoms and other test results.

LIMITATIONS

- As cultures containing AFB and so forth are employed for the tests using this product, the operation described in this package insert must be conducted in a biological safety cabinet, and used test plates must be handled with appropriate precautions as they may cause infections.
- The test plate should be used immediately after opening the packaging. When it absorbs moisture, the quality deteriorates and an accurate result cannot be obtained.
- Any clinical diagnosis must be made by an attending physician, in combination with the clinical symptoms and other test results.
- Please use this product following the operational method and precautions described in this package insert. We cannot guarantee results obtained from any other operations and for any other purposes that are not described in the package insert.

PERFORMANCE CHARACTERISTICS

1. Correlation data

1) WHO Strains of *M. tuberculosis*

		Control product		
		Positive	Negative	Total
This product	Positive	70	0	70
	Negative	0	0	0
	Total	70	0	70

Sensitivity : 100% (70/70)
Specificity : – (0/0)
Accuracy : 100% (70/70)

As the above data show, positive results were also obtained with this product for all of the 70 WHO strains of *M. tuberculosis* for which positive results were obtained with the control product.

2) Clinical Isolate of *M. tuberculosis*

		Control product		
		Positive	Negative	Total
This product	Positive	46	0	46
	Negative	0	5 ^{Note}	5
	Total	46	5	51

Sensitivity : 100% (46/46)
Specificity : 100% (5/5)
Accuracy : 100% (51/51)

As the above data show, positive results were also obtained with this product for all of the 46 clinical isolates of *M. tuberculosis* for which positive results were obtained with the control product.

^{Note} A gene analysis for 5 clinical isolates of *M. tuberculosis* for which negative results were obtained by the control product, as well as by this product, revealed that there was mutation in the base sequence of the MPB64 gene, so that those isolates were mutants in which the expression of the MPB64 protein was incomplete.

2. Sensitivity (Detection limit)

The minimum detection limit is 1.2×10^6 CFU/mL.

3. Reactivity

1) Strains demonstrating the production of MPB64

A reaction was observed with the following 8 strains.

<i>M. tuberculosis</i> H37Rv	<i>M. tuberculosis</i> H37Ra	<i>M. africanum</i>
<i>M. bovis</i> deer	<i>M. microti</i>	<i>M. bovis</i> BCG-Tokyo
<i>M. bovis</i> BCG-Russia	<i>M. bovis</i> BCG-Moreau	

2) Strains demonstrating no production of MPB64

No reaction was observed with the following 2 strains.

<i>M. bovis</i> BCG-Pasteur	<i>M. bovis</i> BCG-Tice
-----------------------------	--------------------------

4. Cross reactivity

The reactivity with the following non-tuberculous mycobacteria was checked and an absence of cross reactivity has been confirmed with all of them.

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INTERFERING SUBSTANCES

Clinical tests have been carried out for this product using cultures from sputum, bronchial lavage fluid, pleural effusion, gastric juice and purulent fluid as samples, and no interference in the test results by clinical samples has so far been observed. In addition, we have used the following media for AFB, and no interference in the test results by the media was observed:

- Egg-based media: 3% Ogawa media, 2% Ogawa media, 1% Ogawa media, Lowenstein-Jensen (LJ) media
- Agar media: Middle Brook 7H10 agar media, Middle Brook 7H11 agar media
- Liquid media: Middle Brook 7H9 liquid media, Dubos liquid media, Kirchner media, Sauton's media

As explained above, no impact from different kinds of clinical samples or culture media has so far been found in the study. However, it is not known whether or not any other substances present in samples may affect test results.

REFERENCES

- 1) Mori T. Current status of tuberculosis [in Japanese]. *J Clin Lab Med.* 1999;43:491-498.
- 2) Shimouchi A. Tuberculosis in overseas [in Japanese]. *Clin Microbiol.* 1997;24:9-14.
- 3) Tuberculosis and Infectious Diseases Control Division, Health Service Bureau, Ministry of Health, Labour and Welfare. *Statistics of Tuberculosis in 2004* [in Japanese]. Japan Anti-Tuberculosis Association;2004.
- 4) Koyama A, et al. *Non-tuberculous Mycobacterium Infection* [in Japanese]. Japan Anti-Tuberculosis Association;1998. JATA Books No.11.

- 5) Okada J. Acid fast bacillus [in Japanese]. *J Clin Lab Med.* 1996;40:589-594.
- 6) Kanai M, et al. Clinical bacterial test. In: Kanai M, et al. *Kanai's Manual of Clinical Laboratory Medicine*, 30th ed [in Japanese]. Kanehara & Co;1993.
- 7) Saito H, Abe C, et al. *Acid Fast Bacillus Tests* [in Japanese]. Ishiyaku Publishers;1997.
- 8) Nagasawa M, et al. Identification of the *M. tuberculosis* complex and *M. avium* complex by chemiluminescence labeled DNA probe [in Japanese]. *J Clin Lab Med.* 1992;36:197-200.
- 9) Kusaba K, et al. Nucleic acid identification from culture isolation [in Japanese]. *J Clin Lab Med.* 1999;43:527-531.
- 10) Goto M. Nucleic acid identification from direct clinical specimen [in Japanese]. *J Clin Lab Med.* 1999;43:533-539.
- 11) Abe C. Application of nucleic acid amplification to detection of tuberculosis bacteria and diagnosis [in Japanese]. *Exp Med.* 1887;15:805-807.
- 12) Harboe M, et al. Properties of proteins MPB64, MPB70, and MPB80 of *Mycobacterium bovis* BCG. *Infect Immun.* 1986;52:293-302.
- 13) Oettinger T, et al. Cloning and B-cell-epitope mapping of MPT64 from *Mycobacterium tuberculosis* H37Rv. *Infect Immun.* 1994;62:2058-2064.
- 14) Yamaguchi R, et al. Cloning and characterization of the gene for immunogenic protein MPB64 of *Mycobacterium bovis* BCG. *Infect Immun.* 1986;57:283-288.
- 15) Andersen AB, et al. MPB64 possesses "tuberculosis-complex"-specific B- and T-cell epitopes. *Scand J Immunol.* 1991;34:365-372.
- 16) Li H, et al. Evidence for absence of the MPB64 gene in some substrains of *Mycobacterium bovis* BCG. *Infect Immun.* 1993;61:1730-1734.
- 17) Abe C, et al. Simple and rapid identification of the *Mycobacterium tuberculosis* complex by immunochromatographic assay using anti-MPB64 monoclonal antibodies. *J Clin Microbiol.* 1999;37:3693-3697.
- 18) Onishi H, et al. Rapid identification of tuberculosis bacteria by the anti-MPB 64 antibody based immunochromatography [in Japanese]. *J Jpn Soc Clin Microbiol.* 1999;9:228-233.
- 19) Furuhashi Y, et al. Rapid identification of the tuberculosis complex by BACTEC MGIT 960 and MPB64 immunochromatography [in Japanese]. In: Abstracts from the 12th General Meeting for the Association for the Rapid Method and Automation in Microbiology 1999;53.
- 20) Hasegawa N, et al. New simple and rapid test for culture confirmation of *Mycobacterium tuberculosis* complex: a multicenter study. *J Clin Microbiol.* 2002;40:908-912.
- 21) Hasegawa M, et al. Evaluation of rapid identification method for *Mycobacterium tuberculosis* complex using the immunochromatographic slide test kit [in Japanese]. *J Jpn Assoc Infect Dis.* 2003;77:110-115.
- 22) Hirano K, et al. Mutation including IS6110 insertion in the gene encoding the MPB64 protein of Capilia TB negative *Mycobacterium tuberculosis* isolates. *J Clin Microbiol.* 2004;42:390-392.

INQUIRES

FAX : +81-558-76-0022

ECREP Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands

GLOSSARY OF SYMBOLS

	CE Marking (European directive 98/79/EC on <i>in vitro</i> diagnostic medical devices)		Authorized representative in the European Community
	<i>In vitro</i> diagnostic medical device		Do not reuse
	Temperature limitation		Manufacturer/Manufactured by
	Use by YYYY-MM		Consult instructions for use
	Batch code		Caution, consult accompanying documents.
	Catalog number		Keep away from sunlight
	Contents sufficient for <n> tests		Fragile, handle with care
	Open here		

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REFERENZEN

- 1) Mori T. Current status of tuberculosis [in Japanese]. *J Clin Lab Med.* 1999;43:491-498.
- 2) Shimouchi A. Tuberculosis overseas [in Japanese]. *Clin Microbiol.* 1997;24:9-14.
- 3) Tuberculosis and Infectious Diseases Control Division, Health Service Bureau, Ministry of Health, Labour and Welfare. *Statistics of Tuberculosis in 2004* [in Japanese]. Japan Anti-Tuberculosis Association; 2004.
- 4) Koyama A, et al. *Non-tuberculous Mycobacterium Infection* [in Japanese]. Japan Anti-Tuberculosis Association; 1998. JATA Books No.11.
- 5) Okada J. Acid fast bacillus [in Japanese]. *J Clin Lab Med.* 1996;40:589-594.
- 6) Kanai M, et al. Clinical bacterial test. In: Kanai M, et al. *Kanai's Manual of Clinical Laboratory Medicine*, 30th ed [in Japanese]. Kanehara & Co; 1993.
- 7) Saito H, Abe C, et al. *Acid Fast Bacillus Tests* [in Japanese]. Ishiyaku Publishers; 1997.
- 8) Nagasawa M, et al. Identification of the *M. tuberculosis* complex and *M. avium* complex by chemiluminescence labeled DNA probe [in Japanese]. *J Clin Lab Med.* 1992;36:197-200.
- 9) Kusaba K, et al. Nucleic acid identification from culture isolation [in Japanese]. *J Clin Lab Med.* 1999;43:527-531.
- 10) Goto M. Nucleic acid identification from direct clinical specimen [in Japanese]. *J Clin Lab Med.* 1999;43:533-539.
- 11) Abe C. Application of nucleic acid amplification to detection of tuberculosis bacteria and diagnosis [in Japanese]. *Exp Med.* 1887;15:805-807.
- 12) Harboe M, et al. Properties of proteins MPB64, MPB70, and MPB80 of *Mycobacterium bovis* BCG. *Infect Immun.* 1986;52:293-302.
- 13) Oettinger T, et al. Cloning and B-cell-epitope mapping of MPT64 from *Mycobacterium tuberculosis* H37Rv. *Infect Immun.* 1994;62:2058-2064.
- 14) Yamaguchi R, et al. Cloning and characterization of the gene for immunogenic protein MPB64 of *Mycobacterium bovis* BCG. *Infect Immun.* 1986;57:283-288.
- 15) Andersen AB, et al. MPB64 possesses "tuberculosis-complex"-specific B- and T-cell epitopes. *Scand J Immunol.* 1991;34:365-372.
- 16) Li H, et al. Evidence for absence of the MPB64 gene in some substrains of *Mycobacterium bovis* BCG. *Infect Immun.* 1993;61:1730-1734.
- 17) Abe C, et al. Simple and rapid identification of the *Mycobacterium tuberculosis* complex by immunochromatographic assay using anti-MPB64 monoclonal antibodies. *J Clin Microbiol.* 1999;37:3693-3697.
- 18) Onishi H, et al. Rapid identification of tuberculosis bacteria by the anti- MPB 64 antibody based immunochromatography [in Japanese]. *J Jpn Soc Clin Microbiol.* 1999;9:228-233.
- 19) Furuhashi Y, et al. Rapid identification of the tuberculosis complex by BACTEC MGIT 960 and MPB64 immunochromatography [in Japanese]. In: Abstracts from the 12th General Meeting for the Association for the Rapid Method and Automation in Microbiology 1999;53.
- 20) Hasegawa N, et al. New simple and rapid test for culture confirmation of *Mycobacterium tuberculosis* complex: a multicenter study. *J Clin Microbiol.* 2002;40:908-912.
- 21) Hasegawa M, et al. Evaluation of rapid identification method for *Mycobacterium tuberculosis* complex using the immunochromatographic slide test kit [in Japanese]. *J Jpn Assoc Infect Dis.* 2003;77:110-115.
- 22) Hirano K, et al. Mutation including IS6110 insertion in the gene encoding the MPB64 protein of Capilia TB negative *Mycobacterium tuberculosis* isolates. *J Clin Microbiol.* 2004;42:390-392.

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REFERENCES

- 1) Mori T. Current status of tuberculosis [in Japanese]. *J Clin Lab Med.* 1999;43:491-498.
- 2) Shimouchi A. Tuberculosis in overseas [in Japanese]. *Clin Microbiol.* 1997;24:9-14.
- 3) Tuberculosis and Infectious Diseases Control Division, Health Service Bureau, Ministry of Health, Labour and Welfare. *Statistics of Tuberculosis in 2004* [in Japanese]. Japan Anti-Tuberculosis Association;2004.
- 4) Koyama A, et al. *Non-tuberculous Mycobacterium Infection* [in Japanese]. Japan Anti-Tuberculosis Association;1998. JATA Books No.11.
- 5) Okada J. Acid fast bacillus [in Japanese]. *J Clin Lab Med.* 1996;40:589-594.
- 6) Kanai M, et al. Clinical bacterial test. In: Kanai M, et al. *Kanai's Manual of Clinical Laboratory Medicine*, 30th ed [in Japanese]. Kanehara & Co;1993.
- 7) Saito H, Abe C, et al. *Acid Fast Bacillus Tests* [in Japanese]. Ishiyaku Publishers;1997.
- 8) Nagasawa M, et al. Identification of the *M. tuberculosis* complex and *M. avium* complex by chemiluminescence labeled DNA probe [in Japanese]. *J Clin Lab Med.* 1992;36:197-200.
- 9) Kusaba K, et al. Nucleic acid identification from culture isolation [in Japanese]. *J Clin Lab Med.* 1999;43:527-531.
- 10) Goto M. Nucleic acid identification from direct clinical specimen [in Japanese]. *J Clin Lab Med.* 1999;43:533-539.
- 11) Abe C. Application of nucleic acid amplification to detection of tuberculosis bacteria and diagnosis [in Japanese]. *Exp Med.* 1887;15:805-807.
- 12) Harboe M, et al. Properties of proteins MPB64, MPB70, and MPB80 of *Mycobacterium bovis* BCG. *Infect Immun.* 1986;52:293-302.
- 13) Oettinger T, et al. Cloning and B-cell-epitope mapping of MPT64 from *Mycobacterium tuberculosis* H37Rv. *Infect Immun.* 1994;62:2058-2064.
- 14) Yamaguchi R, et al. Cloning and characterization of the gene for immunogenic protein MPB64 of *Mycobacterium bovis* BCG. *Infect Immun.* 1986;57:283-288.
- 15) Andersen AB, et al. MPB64 possesses "tuberculosis-complex"-specific B-and T-cell epitopes. *Scand J Immunol.* 1991;34:365-372.
- 16) Li H, et al. Evidence for absence of the MPB64 gene in some substrains of *Mycobacterium bovis* BCG. *Infect Immun.* 1993;61:1730-1734.
- 17) Abe C, et al. Simple and rapid identification of the *Mycobacterium tuberculosis* complex by immunochromatographic assay using anti-MPB64 monoclonal antibodies. *J Clin Microbiol.* 1999;37:3693-3697.
- 18) Onishi H, et al. Rapid identification of tuberculosis bacteria by the anti-MPB 64 antibody based immunochromatography [in Japanese]. *J Jpn Soc Clin Microbiol.* 1999;9:228-233.
- 19) Furuhashi Y, et al. Rapid identification of the tuberculosis complex by BACTEC MGIT 960 and MPB64 immunochromatography [in Japanese]. In: Abstracts from the 12th General Meeting for the Association for the Rapid Method and Automation in Microbiology 1999;53.
- 20) Hasegawa N, et al. New simple and rapid test for culture confirmation of *Mycobacterium tuberculosis* complex: a multicenter study. *J Clin Microbiol.* 2002;40:908-912.
- 21) Hasegawa M, et al. Evaluation of rapid identification method for *Mycobacterium tuberculosis* complex using the immunochromatographic slide test kit [in Japanese]. *J Jpn Assoc Infect Dis.* 2003;77:110-115.
- 22) Hirano K, et al. Mutation including IS6110 insertion in the gene encoding the MPB64 protein of Capilia TB negative *Mycobacterium tuberculosis* isolates. *J Clin Microbiol.* 2004;42:390-392.

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REFERENCES

- 1) Mori T. Current status of tuberculosis [in Japanese]. *J Clin Lab Med.* 1999; 43:491-498.
- 2) Shimouchi A. Tuberculosis in overseas [in Japanese]. *Clin Microbiol.* 1997; 24:9-14.
- 3) Tuberculosis and Infectious Diseases Control Division, Health Service Bureau, Ministry of Health, Labour and Welfare. *Statistics of Tuberculosis in 2004* [in Japanese]. Japan Anti-Tuberculosis Association;2004.
- 4) Koyama A, et al. *Non-tuberculous Mycobacterium Infection* [in Japanese] Japan Anti-Tuberculosis Association;1998. JATA Books No.11.
- 5) Okada J. Acid fast bacillus [in Japanese]. *J Clin Lab Med.* 1996;40:589-594.
- 6) Kanai M, et al. Clinical bacterial test. In: Kanai M, et al. *Kanai's Manual of Clinical Laboratory Medicine*, 30th ed [in Japanese]. Kanehara & Co;1993.
- 7) Saito H, Abe C, et al. Acid Fast Bacillus Tests [in Japanese]. Ishiyaku Publishers;1997.
- 8) Nagasawa M, et al. Identification of the *M. tuberculosis* complex and *M. avium* complex by chemiluminescence labeled DNA probe [in Japanese].*J Clin Lab Med.* 1992;36:197-200.
- 9) Kusaba K, et al. Nucleic acid identification from culture isolation [in Japanese]. *J Clin Lab Med.* 1999;43:527-531
- 10) Goto M. Nucleic acid identification from direct clinical specimen [in Japanese]. *J Clin Lab Med.* 1999;43:533-539.
- 11) Abe C.Application of nucleic acid amplification to detection of tuberculosis bacteria and diagnosis [in Japanese]. *Exp Med.* 1887;15:805-807.
- 12) Harboe M, et al. Properties of proteins MPB64, MPB70, and MPB80 of *Mycobacterium bovis* BCG. *Infect Immun.* 1986;52:293-302.
- 13) Oettinger T,et al.Cloning and B-cell-epitope mapping of MPT64 from *Mycobacterium tuberculosis* H37Rv. *Infect Immun.* 1994;62:2058-2064.
- 14) Yamaguchi R,et al.Cloning and characterization of the gene for immunogenic protein MPB64 of *Mycobacterium bovis* BCG. *Infect Immun.* 1986;57:283- 288.
- 15) Andersen AB, et al. MPB64 possesses "tuberculosis-complex"-specific B-and T-cell epitopes. *Scand J Immunol.* 1991;34:365-372.
- 16) Li H, et al. Evidence for absence of the MPB64 gene in some substrains of *Mycobacterium bovis* BCG. *Infect Immun.* 1993;61:1730-1734.
- 17) Abe C, et al. Simple and rapid identification of the *Mycobacterium tuberculosis* complex by immunochromatographic assay using anti-MPB64 monoclonal antibodies. *J Clin Microbiol.* 1999;37:3693-3697.
- 18) Onishi H, et al. Rapid identification of tuberculosis bacteria by the anti-MPB 64 antibody based immunochromatography [in Japanese]. *J Jpn Soc Clin Microbiol.* 1999;9:228-233.
- 19) Furuhashi Y, et al. Rapid identification of the tuberculosis complex by BACTEC MGIT 960 and MPB64 immunochromatography [in Japanese]. In: Abstracts from the 12th General Meeting for the Association for the Rapid Method and Automation in Microbiology 1999;53.
- 20) Hasegawa N, et al. New simple and rapid test for culture confirmation of *Mycobacterium tuberculosis* complex: a multicenter study. *J Clin Microbiol.* 2002;40:908-912.
- 21) Hasegaw M, et al. Evaluation of rapid identification method for *Mycobacterium tuberculosis* complex using the immunochromatographic slide test kit [in Japanese]. *J Jpn Assoc Infect Dis.* 2003;77:110-115.
- 22) Hirano K, et al. Mutation including IS6110 insertion in the gene encoding the MPB64 protein of Capilia TB negative *Mycobacterium tuberculosis* isolates. *J Clin Microbiol.* 2004;42:390-392.

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REFERENCES

- 1) Mori T. Current status of tuberculosis [in Japanese]. *J Clin Lab Med.* 1999; 43:491-498.
- 2) Shimouchi A. Tuberculosis in overseas [in Japanese]. *Clin Microbiol.* 1997; 24:9-14.
- 3) Tuberculosis and Infectious Diseases Control Division, Health Service Bureau, Ministry of Health, Labour and Welfare. *Statistics of Tuberculosis in 2004* [in Japanese]. Japan Anti-Tuberculosis Association;2004.
- 4) Koyama A, et al. *Non-tuberculous Mycobacterium Infection* [in Japanese] Japan Anti-Tuberculosis Association;1998. JATA Books No.11.
- 5) Okada J. Acid fast bacillus [in Japanese]. *J Clin Lab Med.* 1996;40:589-594.
- 6) Kanai M, et al. Clinical bacterial test. In: Kanai M, et al. *Kanai's Manual of Clinical Laboratory Medicine*, 30th ed [in Japanese]. Kanehara & Co;1993.
- 7) Saito H, Abe C, et al. Acid Fast Bacillus Tests [in Japanese]. Ishiyaku Publishers;1997.
- 8) Nagasawa M, et al. Identification of the *M. tuberculosis complex* and *M. avium complex* by chemiluminescence labeled DNA probe [in Japanese]. *J Clin Lab Med.* 1992;36:197-200.
- 9) Kusaba K, et al. Nucleic acid identification from culture isolation [in Japanese]. *J Clin Lab Med.* 1999;43:527-531
- 10) Goto M. Nucleic acid identification from direct clinical specimen [in Japanese]. *J Clin Lab Med.* 1999;43:533-539.
- 11) Abe C. Application of nucleic acid amplification to detection of tuberculosis bacteria and diagnosis [in Japanese]. *Exp Med.* 1987;15:805-807.
- 12) Harboe M, et al. Properties of proteins MPB64, MPB70, and MPB80 of *Mycobacterium bovis* BCG. *Infect Immun.* 1986;52:293-302.
- 13) Oettinger T, et al. Cloning and B-cell-epitope mapping of MPT64 from *Mycobacterium tuberculosis* H37Rv. *Infect Immun.* 1994;62:2058-2064.
- 14) Yamaguchi R, et al. Cloning and characterization of the gene for immunogenic protein MPB64 of *Mycobacterium bovis* BCG. *Infect Immun.* 1986;57:283-288.
- 15) Andersen AB, et al. MPB64 possesses "tuberculosis-complex"-specific B- and T-cell epitopes. *Scand J Immunol.* 1991;34:365-372.
- 16) Li H, et al. Evidence for absence of the MPB64 gene in some substrains of *Mycobacterium bovis* BCG. *Infect Immun.* 1993;61:1730-1734.
- 17) Abe C, et al. Simple and rapid identification of the *Mycobacterium tuberculosis complex* by immunochromatographic assay using anti-MPB64 monoclonal antibodies. *J Clin Microbiol.* 1999;37:3693-3697.
- 18) Onishi H, et al. Rapid identification of tuberculosis bacteria by the anti- MPB 64 antibody based immunochromatography [in Japanese]. *J Jpn Soc Clin Microbiol.* 1999;9:228-233.
- 19) Furuhashi Y, et al. Rapid identification of the tuberculosis complex by BACTEC MGIT 960 and MPB64 immunochromatography [in Japanese]. In: Abstracts from the 12th General Meeting for the Association for the Rapid Method and Automation in Microbiology 1999;53.
- 20) Hasegawa N, et al. New simple and rapid test for culture confirmation of *Mycobacterium tuberculosis complex*: a multicenter study. *J Clin Microbiol.* 2002;40:908-912.
- 21) Hasegawa M, et al. Evaluation of rapid identification method for *Mycobacterium tuberculosis complex* using the immunochromatographic slide test kit [in Japanese]. *J Jpn Assoc Infect Dis.* 2003;77:110-115.
- 22) Hirano K, et al. Mutation including IS6110 insertion in the gene encoding the MPB64 protein of Capilia TB negative *Mycobacterium tuberculosis isolates*. *J Clin Microbiol.* 2004;42:390-392.

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